



*AAPM Computed Tomography Radiation  
Dose Education Slides  
Philips Version*

Many of the terms used in these slides can  
be found in the CT Terminology Lexicon

[http://www.aapm.org/pubs/CTProtocols/  
documents/CTTerminologyLexicon.pdf](http://www.aapm.org/pubs/CTProtocols/documents/CTTerminologyLexicon.pdf)

Last updated: 18 November 2013



## *Disclaimer*

---

The individuals involved in the production of this document along with the AAPM, ACR, ASRT, MITA, the FDA and the Image Gently Alliance and the Image Wisely Alliance do not endorse any of the specific products mentioned in these slides. Any screenshots of particular products are included for instructional purposes only.



## Disclaimer

---

- Screen captures are **examples** of a common (or latest) software version only and all software versions are not represented
- The information contained herein is current as of the date shown on the title slide
- The master version of these slides is located at:
  - <http://www.aapm.org/pubs/CTProtocols/documents/EducationSlides.pptx>
- Modification of the content of these slides is **not allowed**.



## *Vendor Specific Slide Details*

---

- The presence of a vendor name in the title of the slide indicates that the slide is vendor specific slide
- White text is used throughout to indicate vendor specific language
- An example of a vendor specific slide follows



## Vendor: *Generic Parameter/Topic Name*

---

Vendor Specific Name

Vendor screen  
capture of how the  
acquisition  
parameter is set  
or how  
information on the  
topic is displayed

Text describing acquisition  
parameter or topic



## Motivation

- These slides are provided to aid in understanding the factors that affect radiation dose in CT studies
- Image patients **wisely** and **gently**
  - A CT study should use as little radiation as possible, while still meeting the image quality needs of the exam
  - A CT study that is non-diagnostic because the radiation dose is too low may require rescanning the patient – increasing the total patient dose



[imagegently.org](http://imagegently.org)

[imagewisely.org](http://imagewisely.org)





## Outline

---

- What is Dose?
- Acquisition Parameter Settings
- Dose Modulation and Reduction
- Dose Display



## What Is Dose?

- Volume Computed Tomography Dose Index ( $CTDI_{vol}$ ) is a standardized parameter to measure **Scanner Radiation Output**
  - $CTDI_{vol}$  is NOT patient dose
  - $CTDI_{vol}$  is reported in units of mGy for either a 16-cm (for head exams) or 32-cm (for body exams) diameter acrylic phantom
  - For the same technique settings, the  $CTDI_{vol}$  reported for the 16-cm phantom is about twice that of the 32-cm phantom
  - The reported  $CTDI_{vol}$  is based on measurements made by the manufacturer in a factory setting
- In these slides, the term "patient dose" is used to describe the absorbed dose to a patient, while the generic term "dose" refers to  $CTDI_{vol}$

1. Bauhs, J. A., Vrieze, T. J., Primak, A. N., Bruesewitz, M. R., & McCollough, C. H. (2008). CT Dosimetry: Comparison of Measurement Techniques and Devices1. *Radiographics*, 28(1), 245-253. doi:10.1148/rg.281075024
2. McCollough, C. H., Primak, A. N., Braun, N., Kofler, J., Yu, L., & Christner, J. (2009). Strategies for reducing radiation dose in CT. *Radiologic clinics of North America*, 47(1), 27-40.
3. International Electrotechnical Commission. *Medical Electrical Equipment. Part 2–44: Particular requirements for the safety of x-ray equipment for computed tomography*. 2.1. International Electrotechnical Commission (IEC) Central Office; Geneva, Switzerland: 2002. IEC publication No. 60601–2–44.





## *How is $CTDI_{vol}$ related to patient dose?*

---

- $CTDI_{vol}$  is not patient dose
- The relationship between the two depends on many factors, including patient size and composition
- [AAPM Report 204](#) introduces a parameter known as the Size Specific Dose Estimate (SSDE) to allow estimation of patient dose based on  $CTDI_{vol}$  and patient size
- For the same  $CTDI_{vol}$ , a smaller patient will tend to have a higher patient dose than a larger patient

What is Dose?

[http://www.aapm.org/pubs/reports/RPT\\_204.pdf](http://www.aapm.org/pubs/reports/RPT_204.pdf)



## How is $CTDI_{vol}$ related to patient dose?

120 kVp at 200 mAs



32 cm  
Phantom

$CTDI_{vol} = 20 \text{ mGy}$

120 kVp at 200 mAs



32 cm  
Phantom

$CTDI_{vol} = 20 \text{ mGy}$

**Both patients scanned with the same  $CTDI_{vol}$**   
**Patient dose will be higher for the smaller patient**

What is Dose?



## How is $CTDI_{vol}$ related to patient dose?

120 kVp at 100 mAs



32 cm  
Phantom

$CTDI_{vol} = 10 \text{ mGy}$

120 kVp at 200 mAs



32 cm  
Phantom

$CTDI_{vol} = 20 \text{ mGy}$

**Smaller patient scanned with a lower  $CTDI_{vol}$   
Patient doses will be approximately equal**

What is Dose?



## Size Specific Dose Estimate (SSDE)


- AAPM report 204 describes a method to calculate SSDE using  $CTDI_{vol}$
- Conversion factors based on patient size (e.g., AP or lateral width, effective diameter) are provided to **estimate** patient dose for a patient of that size
- However, SSDE is still not the exact patient dose, as factors such as scan length and patient composition may differ from the assumptions used to calculate SSDE
- SSDE is not dose to any specific organ, but rather the mean dose in the center of the scanned volume

What is Dose?




## How is $CTDI_{vol}$ related to patient dose?

120 kVp at 100 mAs



$CTDI_{vol} = 10 \text{ mGy}$   
 $SSDE = 13.2 \text{ mGy}$

120 kVp at 200 mAs



$CTDI_{vol} = 20 \text{ mGy}$   
 $SSDE = 13.2 \text{ mGy}$

**Patients have equivalent SSDE**

What is Dose?



## Why Use $CTDI_{vol}$ ?

- $CTDI_{vol}$  provides information about the amount of radiation used to perform the study
- $CTDI_{vol}$  is a useful index to track across patients and protocols for quality assurance purposes
- $CTDI_{vol}$  can be used as a metric to compare protocols across different practices and scanners when related variables, such as resultant image quality, are also taken in account
- The ACR Dose Index Registry (DIR) allows comparison across institutions of  $CTDI_{vol}$  for similar exam types (e.g., routine head exam)

What is Dose?

1. McCollough, C. H., Leng, S., Yu, L., Cody, D. D., Boone, J. M., & McNitt-Gray, M. F. (2011). CT Dose Index and Patient Dose: They are Not the Same Thing, EDITORIAL, *Radiology* 259(2), 311-316.



## *Dose Length Product*

---

- The Dose Length Product (DLP) is also calculated by the scanner
- DLP is the product of the length of the irradiated scan volume and the average  $CTDI_{vol}$  over that distance
- DLP has units of  $mGy \cdot cm$

What is Dose?



## Useful Concepts/Terms

- The relationships between acquisition parameters and  $CTDI_{vol}$  described in the following slides assume all other parameters are held constant
- The relationship between a parameter and  $CTDI_{vol}$  is often described as **proportional** in some way
  - The symbol  $\propto$  is used to indicate “proportional to”
- Directly proportional means that a change in the parameter results in the same change in  $CTDI_{vol}$ 
  - Example: Doubling the rotation time from 0.5 to 1.0 seconds will double the  $CTDI_{vol}$
- Inversely proportional means that a change in a parameter has the opposite effect on  $CTDI_{vol}$ 
  - Example: Doubling the pitch from 1 to 2 will reduce the  $CTDI_{vol}$  by half





## *Acquisition Parameter Settings*

---

- Acquisition Parameters define the technique that will be used and how the scan will proceed
- Acquisition Parameters are set in the user interface where scans are prescribed
- Changing a single Acquisition Parameter while holding everything else constant will typically affect the  $CTDI_{vol}$  for that scan
- The following slides describe what that affect is for each parameter



## Scan Mode

---

- CT Scanners offer a variety of **Scan Modes** which describe how the table moves during an exam
- **Scan Modes** include
  - Axial
  - Helical or Spiral
  - Dynamic

**The Acquisition Parameters that affect CTDIvol may change amongst different Scan Modes**

Acquisition Parameter Settings



## Dynamic Scan Mode Notes

- In the Dynamic Scan Mode multiple acquisitions covering the same body region are acquired. Examples of these study types include:
  - Perfusion Studies
  - Bolus Tracking Studies
  - Test Bolus Studies
- Dynamic Scans often have large  $CTDI_{vol}$  values because the scanner reports the sum of the  $CTDI_{vol}$  values from each rotation
- The reported  $CTDI_{vol}$  is NOT skin dose or organ dose


Acquisition Parameter Settings

1. Bauhs, J. A., Vrieze, T. J., Primak, A. N., Bruesewitz, M. R., & McCollough, C. H. (2008). CT Dosimetry : Comparison of Measurement Techniques and Devices. *Radiographics*, 28(1), 245-254.
2. Zhang, D., Cagnon, C. H., Villablanca, J. P., McCollough, C. H., Cody, D. D., Stevens, D. M., Zankl, M., et al. (2012). Peak Skin and Eye Lens Radiation Dose From Brain Perfusion CT Based on Monte Carlo Simulation. *American Journal of Roentgenology*, 198(2), 412-417.



# Philips: Scan Mode

## Scan Mode



- Scan mode is defined during Exam card creation
- Dynamic scans, such as Brain Perfusion use '0' scan increment.
- In axial scanning, the scan increment will usually equal the collimation width ( $64 \times 0.625 = 40$  mm), unless the user chooses a larger increment, creating 'gaps' between the axial acquisitions.

Acquisition Parameter Settings



## Philips: Scan Mode

### Scan Mode

6 Brain, Helical

General

Label:

Collimation:  (64x0.625)

Sampled Collimation:

Pitch:  ←

Rotation Time:

Automatic Scan Time

Scan Time:  sec

Auto Voice

Edit before final Recon

- Helical Scans show 'pitch' instead of 'increment' in the 'Show All' acquisition parameters screen.

Acquisition Parameter Settings



## *Table Feed/Increment*

---

- Is the movement of the table through the bore of the scanner over a full 360 degree rotation
- Units: millimeters/rotation or millimeters/second
- The parameter is known both as **Table Feed** (helical/spiral acquisition) & **Table Increment** (axial acquisition)

**Table Feed may affect  $CTDI_{vol}$  through its inclusion in Pitch (discussed later)**

Acquisition Parameter Settings



# Philips: Table Feed/Increment

## Table Speed

Study ID	49	Slice Increment	-3.00 mm
Content Date	3 Dec, 2012	Pitch	0.67
Content Time	9:56:22.16	Rotation Time	1.00 Sec
Patient Position	Head First Supine	Field of View	302.0 mm
Scan Options	HELIX	Voltage	120 kV
Slice Location	108.2 mm	Current	435 mA
Table Height	85	mAs	649 mAs/Slice
Table Speed	26.8 mm/sec	Dose Modulation	3D MODULATION
Scan Length	130.1 mm	SP Filter	Yes
Scan Time	7.20 Sec	Adaptive Filter	YES
Image Type	AXIAL	Filter	B
		Image Matrix	512 x 512
		Center X, Y	0.0, 0.0

- Increment for Helical scans represents image reconstruction increment.
- Table Speed (shown on Image Parameters Page in Viewer) is reported in mm/sec
- Table Speed = Collimation \* Pitch / Rotation Time

Acquisition Parameter Settings



## Detector Configuration

- Is the combination of the number of data channels and the width of the detector associated with each data channel
- The **Detector Configuration** determines the Beam Width or Beam Collimation ( $nT$ ), which is the number of channels ( $n$ ) times the detector width associated with each data channel ( $T$ )
- For a selected detector width per data channel, a smaller total Beam Collimation usually has a higher  $CTDI_{vol}$  than a larger Beam Collimation
  - Example: On a 16 slice scanner with a detector width per channel of 1.25 mm, a collimation of  $4 \times 1.25\text{mm}$  is generally less dose efficient than a collimation of  $16 \times 1.25\text{mm}$

**Users should monitor  $CTDI_{vol}$  values when changing detector configuration**


Acquisition Parameter Settings





# Philips: *Detector Configuration*

## Collimation




- Collimation is selected by clicking on the Collimation dropdown menu on the Show All acquisition parameters screen

Acquisition Parameter Settings



## Detector Configuration



Acquisition Parameter Settings



## *Pitch*

---

- Is the Table Feed per gantry rotation divided by the beam width/collimation
- *Pitch* is the ratio of two distances and therefore has no units
- Users should monitor other parameters when changing *Pitch*. The scanner may or may not automatically compensate for changes in *Pitch* (for example, by changing the tube current) to maintain the planned  $CTDI_{vol}$ .

**$CTDI_{vol} \propto 1/Pitch$ :**

**Hitachi, Toshiba (no AEC)**

**$CTDI_{vol}$  independent of *Pitch*:**

**GE, Siemens, Philips, Neusoft, Toshiba (AEC)**

Acquisition Parameter Settings



## Pitch


---

- $CTDI_{vol}$  may not change in the expected manner if the scanner automatically adjust other parameters when the pitch is changed
- The relationships between  $CTDI_{vol}$  and pitch for the different vendors are described below
  - $CTDI_{vol}$  inversely proportional to change in pitch: Hitachi and NeuroLogica
  - $CTDI_{vol}$  constant when pitch is changed due to changes to other parameters: GE, Neusoft, Philips and Siemens
  - The relationship between  $CTDI_{vol}$  and pitch depends on scan mode or Software version: Toshiba



# Philips: Pitch

## Pitch




- Pitch is selected by clicking on the Pitch dropdown menu on the Show All acquisition parameters screen.
- On Philips scanners, the mA is adjusted to maintain mAs when pitch is adjusted.

Acquisition Parameter Settings



## Pitch



**Pitch < 1**  
Beam Width has some overlap at each view angle from rotation to rotation

**Pitch = 1**  
No overlap of Beam Width at each view angle and no view angles not covered at certain table positions

**Pitch > 1**  
Some view angles are not covered by the beam width at certain table positions

Acquisition Parameter Settings



## *Exposure Time per Rotation*

- Is the length of time, in seconds, that the X-ray beam is “on” during a gantry rotation
  - It takes into account the gantry rotation time and angular acquisition range
- Units: seconds
- Users should monitor other parameters when changing **Exposure Time per Rotation**. The scanner may or may not automatically compensate for changes in **Exposure Time per Rotation**(for example, by changing the tube current)

**$CTDI_{vol} \propto$  Exposure Time per Rotation**

**Hitachi, NeuroLogica, Toshiba (no AEC)**

**$CTDI_{vol}$  independent of Exposure Time per Rotation:**

**GE, Siemens, Philips, Neusoft, Toshiba (AEC)**

Acquisition Parameter Settings



## *Exposure Time per Rotation*

---


- $CTDI_{vol}$  may not change in the expected manner if the scanner automatically adjust other parameters when the exposure time per rotation is changed
- The relationships between  $CTDI_{vol}$  and exposure time per rotation for the different vendors are described below
  - $CTDI_{vol}$  proportional to change in parameter: Hitachi and NeuroLogica
  - $CTDI_{vol}$  constant when the parameter is changed due to changes to other parameters: GE, Neusoft, Philips and Siemens
  - The relationship between  $CTDI_{vol}$  and the parameter depends on scan mode or Software version: Toshiba





## Philips: Exposure Time per Rotation

### Rotation Time



- Rotation Time is selected by clicking on the Rotation Time dropdown menu on the Show All acquisition parameters screen.
- When Automatic Scan time is selected, the system automatically selects the rotation time and pitch based on the DoseRight Index parameters.

Acquisition Parameter Settings



## *Tube Current*

---

- Determines the number of electrons accelerated across the x-ray tube per unit time
- Units: milliAmperes (mA)
- $CTDI_{vol}$  is directly proportional to **Tube Current**

$$CTDI_{vol} \propto \text{Tube Current}$$

Acquisition Parameter Settings



# Philips: Tube Current

## mA

The screenshot displays the Philips CT console interface. On the left, a tree view shows the scan protocol: 1. Survey, Frontal; 2. Chest, Helical (selected). Below this, the 'Scan Parameters' section shows: KV: 120, mAs: 200 (483), Scan Time: 7.4 sec. The main area is titled 'Overview' and 'All Parameters'. It is divided into several sections: 'General' (Label, Collimation: Auto (640.625), Pitch: 1.015, Rotation Time: 0.4, Automatic Scan Time checked, Scan Time: 7.4 sec, Auto Voice: None, Edit before final Recon checked), 'Dose Management' (DoseRight unchecked, KV: 120, mAs: 200 (483), CTDIvol: 13.1 mGy, DLP: 940.4 mGy\*cm, Accumulated CTDI: 13.1 mGy, Accumulated DLP: 940.4 mGy\*cm, Phantom: CTDIvol Body Phantom 32cm), 'Injection' (Injection Parameters..., Contrast, Trigger: Manual), and 'Dose Notification' (Dose Notification Value CTDIvol: 50 mGy, Dose Notification Value DLP: None mGy\*cm). A 'Show All...' button is at the bottom left. Two red arrows point to the 'mAs' field in the 'Scan Parameters' section and the 'mAs' field in the 'Dose Management' section.

- mAs can be viewed or modified in Scan Parameters or the Show All scan parameters (light gray text) and can be modified by changing the pitch, mAs and/or rotation time.



## Tube Potential

- Is the electrical potential applied across the x-ray tube to accelerate electrons toward the target material
- Units: kiloVolts (kV or kVp)
- $\text{CTDI}_{\text{vol}}$  is **approximately** proportional to the square of the percentage change in **Tube Potential**

$$\text{CTDI}_{\text{vol}} \propto \left( \frac{kV_{\text{new}}}{kV_{\text{old}}} \right)^n \quad n \approx 2 \text{ to } 3$$

Acquisition Parameter Settings



# Philips: Tube Potential

## KV

The screenshot shows the Philips CT console interface. On the left, the 'Scan Parameters' section is visible, showing a dropdown menu for 'KV' set to 120. A red arrow points to this dropdown. On the right, the 'Overview' tab is active, showing the 'All Parameters' section. In the 'Dose Management' section, the 'KV' dropdown is also set to 120, with another red arrow pointing to it. The 'Dose Management' section also displays 'mAs: (mA) 200 (483)', 'CTDIvol: 13.1 mGy', and 'DLP: 940.4 mGy\*cm'. The 'Dose Notification' section shows 'Dose Notification Value CTDIvol: 50 mGy' and 'Dose Notification Value DLP: None mGy\*cm'.

- kV can be modified in Scan Parameters or the Show All scan parameters.



## *Tube Current Time Product*

---

- Is the product of Tube Current and the Exposure Time per Rotation
- Units: milliAmpere-seconds (mAs)
- $CTDI_{vol}$  is directly proportional to Tube Current Time Product

$$CTDI_{vol} \propto \text{Tube Current Time Product}$$

Acquisition Parameter Settings



## *Effective Tube Current Time Product*

---

- Is the product of the Tube Current and the Exposure Time per Rotation divided by the Pitch
- Units: milliAmpere-Seconds (mAs)
- $CTDI_{vol}$  is directly proportional to Effective Tube Current Time Product


$$CTDI_{vol} \propto \text{Effective Tube Current Time Product}$$

Acquisition Parameter Settings



# Philips: Effective tube Current Time Product

## mAs



- mAs can be modified in Scan Parameters or the Show All scan parameters, or by raising or lowering the DRI (DoseRight Index ) when DoseRight is enabled.

Acquisition Parameter Settings





## *Field Of Measurement*

---

- Is the diameter of the primary beam in the axial plane at the gantry iso-center
- Units: millimeters (mm)
- $CTDI_{vol}$  may decrease with a decrease in the Field of Measurement
  - The relationship is vendor specific


**Users should monitor the  $CTDI_{vol}$  values when changing the Field of Measurement**

Acquisition Parameter Settings



## Philips: Field of Measurement

### Field of View



- FOV (field of view) can be modified by changing the value in the Exam Card Manager screen when generating protocols.
- Reconstruction FOV can be modified during planning by adjusting the result FOV.

Acquisition Parameter Settings



## *Beam Shaping Filter*

---

- Is the scanner component that modifies the energy spectrum and spatial distribution of the primary beam
- **Beam Shaping** may include a bow tie filter and/or flat filters
- $CTDI_{vol}$  is affected by a change in **Beam Shaping Filters**
  - The relationship is vendor and filter specific

**Users should monitor  $CTDI_{vol}$  values when changing the Beam Shaping Filter**

Acquisition Parameter Settings



## Acquisition Parameter Settings Summary

Parameter	Relationship to $CTDI_{vol}$
Scan Mode	Changes in the Scan Mode may affect $CTDI_{vol}$
Table Feed/Increment	Table Feed affects $CTDI_{vol}$ through its inclusion in Pitch
Detector Configuration	Decreasing the Beam Collimation typically, but not always, increases the $CTDI_{vol}$
Pitch	$CTDI_{vol}$ relationship to pitch is vendor dependent
Exposure Time Per Rotation	$CTDI_{vol}$ relationship to exposure time per rotation is vendor dependent
Tube Current	$CTDI_{vol} \propto$ Tube Current
Tube Potential	$CTDI_{vol} \propto (kVp_1/kVp_2)^n$ $n \sim 2$ to $3$
Tube Current Time Product	$CTDI_{vol} \propto$ Tube Current Time Product
Effective Tube Current Time Product	$CTDI_{vol} \propto$ Effective Tube Current Time Product
Field of Measurement	Changes in the Field of Measurement may affect $CTDI_{vol}$
Beam Shaping Filter	Changes in the Beam Shaping Filter may affect $CTDI_{vol}$



## *Dose Modulation and Reduction*

---

- Many CT scanners automatically adjust the technique parameters (and as a result the  $CTDI_{vol}$ ) to achieve a desired level of image quality and/or to reduce dose
- Dose Modulation and Reduction techniques vary by scanner manufacturer, model and software version



## Automatic Exposure Control (AEC)

- Automatically adapts the Tube Current or Tube Potential according to patient attenuation to achieve a specified image quality
  - Automatic adjustment of Tube Current may not occur when Tube Potential is changed
  - **Centering the patient in the gantry is VITAL for most AEC systems**
- AEC aims to deliver a specified image quality across a range of patient sizes. It tends to increase  $CTDI_{vol}$  for large patients and decrease it for small patients relative to a reference patient size


**The use of Automatic Exposure Control may decrease or increase  $CTDI_{vol}$  depending on the patient size and body area imaged and image quality requested**

Dose Modulation and Reduction



# Philips: Automatic Exposure Control (AEC)

## DoseRight



- DoseRight (AEC) can be modified in Show All scan parameters .
- DoseRight Index (DRI) can be disabled entirely by unchecking the DoseRight box. DRI is an image quality reference parameter.

Dose Modulation and Reduction



## *Image Quality Reference Parameter*

---

- Is the AEC parameter that is set by the user to define the desired level of image quality
- Changing the Image Quality Reference Parameter will affect the  $CTDI_{vol}$

**The effect on  $CTDI_{vol}$  when changing the Image Quality Reference Parameter is vendor dependent**


Dose Modulation and Reduction





# Philips: Image Quality Reference Parameter

## DoseRight Index



DoseRight Index (DRI) is an image quality reference parameter that can be modified in Scan Parameters or Show All scan parameters by increasing or decreasing the DRI value. Dose Modulation and Reduction



## *Philips: Image Quality Reference Parameter*

---

- A change in the Image Quality Reference Parameter will affect the  $CTDI_{vol}$
- Setting the parameter for “increased” image quality (e.g., lower noise) will result in more dose
  - Dose Right Index: Increasing the parameter will result in an increase in the  $CTDI_{vol}$
- Setting the parameter for “decreased” image quality (e.g., more noise) will result in less dose
  - Dose Right Index: Decreasing the parameter will result in a decrease in the  $CTDI_{vol}$

Dose Modulation and Reduction



## *Angular Tube Current Modulation*

---

- Is an AEC feature that adjusts the Tube Current as the x-ray tube rotates around the patient to compensate for attenuation changes with view angle
- **Angular Tube Current Modulation** is used to adjust the Tube Current to attempt to deliver similar dose to the detector at all view angles


**The use of Angular Tube Current Modulation may decrease or increase  $CTDI_{vol}$  depending on the patient size and body area imaged and image quality requested**

Dose Modulation and Reduction



## Philips: Angular Tube Current Modulation


- Angular Tube Current Modulation (3D Modulation) uses information from
  - 1 or 2 view localizer radiograph





# Philips: Angular tube Current Modulation

## 3D Modulation



- 3D modulation can be used together with DoseRight Index or can be turned off at time of exam (This view is from Exam Card manager).
- DRI (AEC) can be disabled entirely by unchecking the DoseRight box in the active exam. In this case, the default mAs will be assumed.

Dose Modulation and Reduction



## *Longitudinal Tube Current Modulation*

---

- Is an AEC feature that adjusts the Tube Current as patient attenuation changes in the longitudinal direction
- The CT Localizer Radiograph is used to estimate patient attenuation


**The use of Longitudinal Tube Current Modulation may decrease or increase  $CTDI_{vol}$  depending on the patient size and body area imaged and image quality requested**

Dose Modulation and Reduction



## Longitudinal Tube Current Modulation

- Longitudinal Tube Current Modulation (Z Modulation) uses information from
  - One view localizer scans




Dose Modulation and Reduction



# Philips: Longitudinal Tube Current Modulation

## Z-Modulation



- Z- modulation can be used together with DoseRight Index or can be turned off at time of exam (This view is from Exam Card manager).
- DRI (AEC) can be disabled entirely by unchecking the DoseRight box in the active exam. In this case, the default mAs will be assumed.

Dose Modulation and Reduction





## *Angular and Longitudinal Tube Current Modulation*


- Is an AEC feature that incorporates the properties of both **Angular and Longitudinal Tube Current Modulation** to
  - Adjust the Tube Current based on the patient's overall attenuation
  - Modulate the Tube Current in the angular (X-Y) and longitudinal (Z) dimensions to adapt to the patient's shape

**The use of Angular and Longitudinal Tube Current Modulation may decrease or increase  $CTDI_{vol}$  depending on the patient size and body area imaged and image quality requested**

Dose Modulation and Reduction



# Angular and Longitudinal Tube Current Modulation




Dose Modulation and Reduction



## Philips: Angular and Longitudinal Tube Current Modulation

### Z-Modulation and 3D Modulation

- 3D Modulation and Z-Modulation are enabled in exam card manager.
- These features can be disabled at the time of exam (this view from Exam Card manager).
- Z-modulation must be enabled when using 3D modulation.



Dose Modulation and Reduction



## *ECG-Based Tube Current Modulation*

---

- Is an AEC feature used with prospectively gated cardiac imaging that adjusts the Tube Current based on the phase within the cardiac cycle
- There are important heart rate considerations to take into account when using prospective gating

**The use of ECG-Based Tube Current Modulation with prospective gating will decrease  $CTDI_{vol}$  compared to retrospective gating**

Dose Modulation and Reduction



## Philips: ECG-Based Tube Current Modulation

### Cardiac DoseRight

**General**

Label:

Scan Type:

Collimation:  (64x0.625)

Pitch:

Rotation Time:

Automatic Scan Time

Scan Time:  sec

Resolution:

Auto Voice

Edit before final Recon

**Dose Management**

DoseRight

DoseRight Index:  Reduce Dose  R D

Z-Modulation

3D Modulation

Cardiac DoseRight

KV:

mAs for Average Adult: (mA)  (371)

Absolute Min mAs:

Absolute Max mAs:

Ref. Size: 79 cm

CTDIvol: 51 mGy DLP: 986 mGy\*cm

Accumulated DLP: 1060.4 mGy\*cm

Phantom: CTDIvol Body Phantom 32cm

**Injection**


**Dose Notification**

- DoseRight Cardiac is enabled by checking the DoseRight Cardiac Box (in exam card manager).
- DRI can also be enabled to take advantage of the dose right index recommendation.

Dose Modulation and Reduction



## ECG-Based Tube Current Modulation





## *Iterative Reconstruction*

---

- Is a feature that uses the information acquired during the scan and repeated reconstruction steps to produce an image with less “noise” or better image quality (e.g., higher spatial resolution or decreased artifacts) than is achievable using standard reconstruction techniques


**The use of Iterative Reconstruction by itself may not decrease  $CTDI_{vol}$ ; with use of Iterative Reconstruction, image quality will change and this may allow a reduction in the  $CTDI_{vol}$  by adjusting the acquisition parameters used for the exam**

Dose Modulation and Reduction



## Philips: Iterative Reconstruction

### iDose<sup>4</sup>




- iDose<sup>4</sup> iterative reconstruction technique can be accessed in the 'Results' Show All parameters. As the iDose<sup>4</sup> level increases, noise reduction increases. The user may generate multiple iDose<sup>4</sup> reconstructions, at various iDose<sup>4</sup> levels, for any given series. Dose Modulation and Reduction





## Philips: Iterative Reconstruction

### IMR



- IMR knowledge based iterative reconstruction technique can be accessed in the 'Results' Show All parameters. As the IMR level increases, noise reduction increases. The user may generate multiple IMR reconstructions, at various IMR levels, and various IMR Image Definitions for any given series.  
Dose Modulation and Reduction



## *Philips: Iterative Reconstruction*

- Iterative Reconstruction using the iDose<sup>4</sup> and IMR iterative reconstruction techniques operates using both the Projection space and Image space data.
- Turning On or Changing the level of iDose<sup>4</sup> or IMR will affect the resultant image quality; it **WILL NOT** affect the  $CTDI_{vol}$  of the scan
- In consultation, the radiologists and medical physicists may adjust the acquisition parameters for studies reconstructed using iDose<sup>4</sup> or IMR based on the imaging task and patient population, dose concerns, and the needs of the interpreting radiologist(s)

Dose Modulation and Reduction



## *Dose Display*

---

- Information about the  $CTDI_{vol}$  planned for each scan is typically displayed before the exam on the user console
- Information about the  $CTDI_{vol}$  delivered by each scan is typically reported in a data page or DICOM structured dose report
- Dose information provided after the exam typically also includes the DLP and the CTDI phantom size. These may also be included in information displayed before the scan.



## *Display of Planned $CTDI_{vol}$*

---

- $CTDI_{vol}$  is displayed before a study is performed based on the selected technique parameters
- It is important to check  $CTDI_{vol}$  before a study is performed to ensure that the output of the scanner is appropriate for the specific patient and diagnostic task

**$CTDI_{vol}$  is displayed for each planned acquisition**

Dose Display



## Philips: Display of Planned $CTDI_{vol}$

The screenshot displays the Philips CT console interface for an 'Abdomen/Pelvis, Helical' scan. The interface is divided into several sections:

- Scan Parameters:** Located in the bottom-left, it shows 'DoseRight Index: 23', 'KV: 120', 'mAs: 163 - 227', and 'Scan Time: 4.4 sec'. Below these, it displays 'CTDIvol: 12.7mGy', 'DLP: 397.6mGy\*cm', and 'Phantom: 32cm'. A 'Show All...' button is present at the bottom.
- General:** Includes 'Collimation: Auto', 'Pitch: 0.891', 'Rotation Time: 0.5', and 'Scan Time: 4.4 sec'. It also has checkboxes for 'Automatic Scan Time', 'Auto Voice', and 'Edit before final Recon'.
- Dose Management:** Features a 'DoseRight' checkbox, 'DoseRight Index' (23), 'KV: 120', 'Average mAs: 193', 'Min mAs: 163', 'Max mAs: 227', and 'Liver Area DoseRight Index: -3'. It also shows 'Ref. Size: 29 cm', 'Patient Size: 23 cm', and 'CTDIvol: 12.7 mGy'.
- Injection:** Includes checkboxes for 'Injection' and 'Contrast', and a 'Trigger: Manual' dropdown.
- Dose Notification:** Shows 'Dose Notification Value CTDIvol: 50 mGy' and 'Dose Notification Value DLP: None mGy\*cm'.

Red arrows point from the text below to the 'CTDIvol: 12.7mGy' and 'CTDIvol: 12.7 mGy' values in the interface.

Planned  $CTDI_{vol}$  is displayed in the Scan Parameters box and Show All scan parameters.

$CTDI_{vol}$  phantom size is displayed.

Dose Display



## *Post Study Data Page*

---

- Following the completion of a study, a **Post Study Data Page** is created that includes information on the delivered  $CTDI_{vol}$  and DLP and the phantom size used to calculate these values
- Information is displayed for each series

Dose Display



# Philips: Post Study Data Page

**CHEST 3.6 mSv**  
 INGENUITY 0/78Y  
 501-2 Dose Information

Philips, Ingenuity CT  
 18 Aug, 2012 10:30:15:35  
 7.1.00

Exam Information						
Study ID	364					
Time	18/08/2012 10:16 AM					
Total DLP	272.6 mGy*cm					

Dose							
#	Scan Label	Scan Mode	mAs	kV	CTDIvol [mGy]	DLP [mGy*cm]	Phantom Type [cm]
1	SURVIEW	Survview	100	0	0.1	0.1	BODY 32 CM
1	SURVIEW	Survview	100	0	0.1	0.1	BODY 32 CM
2	LOCATOR	Station...	30	100	1.5	1.5	BODY 32 CM
3	TRACKER	Station...	30	100	1.5x6	9	BODY 32 CM
4	ACQUISITION	Helical	185	100	7.3	261.9	BODY 32 CM

**PHILIPS**

- This page is displayed at the end of each patient exam and is included as part of the patient record.
- CTDI<sub>vol</sub> is displayed per series along with the DLP. The associated reference phantom type and size is listed for each series.
- A DICOM Dose Structured Report (DICOM SR) output is also generated.

Dose Display



## *Post Study Data Page - $CTDI_{vol}$*

---

- $CTDI_{vol}$  is displayed for each series after a study is performed and is calculated based on the technique factors used to acquire the data
- It is useful to check  $CTDI_{vol}$  after a study is performed to ensure that the output of the scanner was as expected

**$CTDI_{vol}$  is displayed for each completed acquisition**

Dose Display





## *Post Study Data Page - DLP*

---

- DLP is displayed for each series after a study is performed and is calculated based on the technique factors and scan length used

**DLP is displayed for each completed acquisition and is typically summed for all of the acquisitions**

Dose Display



## *Post Study Data Page – CTDI Phantom*

---

- The CTDI Phantom used for each acquisition in the study is typically displayed
- Different phantoms may be used to calculate the  $CTDI_{vol}$  for different acquisitions in the same study (and may vary by vendor)
  - Head and C-Spine Example
    - Body Phantom used to report  $CTDI_{vol}$  for C-Spine portion of exam
    - Head Phantom used to report  $CTDI_{vol}$  for Head portion of exam

Dose Display



## Summing Dose Report Values

- $CTDI_{vol}$  values for separate series are NOT to be summed to give a “total”  $CTDI_{vol}$  for a study
  - This is especially true if the series cover different anatomic regions
- DLP is typically summed over all series in the Post Study Data Page to provide an estimate of the total patient exposure
  - Extreme care should be taken when considering summed DLPs because different phantoms may have been used to calculate the  $CTDI_{vol}$  values used to determine DLP
- A medical physicist should be contacted if patient specific dose estimates are required

Dose Display



## *Dose Notification Levels*

---


- **Notification Levels** may be set on a CT scanner for each series within an exam protocol
- If the planned  $CTDI_{vol}$  is above the **Notification Level** and triggers the notification, the user has the opportunity to edit or confirm the technique settings
- **Notification Levels** may be exceeded when appropriate for a specific patient or diagnostic task (e.g., in very large patients or contrast bolus monitoring scans)

Dose Display



## Philips: Dose Notification Levels

Notification levels are user configurable



- The Dose Notification value is set in the Exam Card (for each series) in Exam Card Manager.

The Dose Notification feature complies with the NEMA XR-25 standard.

Dose Display



## *Dose Alert Levels*

---


- **Dose Alert Levels** require specific action by the operator to continue scanning
- **Dose Alert Levels** are typically much higher than Notification Levels and take into account all series within the exam
- Triggering a **Dose Alert** requires that the operator confirm the protocol and settings are correct by entering in his or her name. Optionally, sites may require that the operator provide a brief explanation in the provided field

Dose Display



## Philips: Dose Alert Levels

Default alert level is 1000 mGy



- The Dose Alert is set in Preferences (Dose Management) and supports Alert levels for both body and head scan types.

The Dose Alert feature complies with the NEMA XR-25 standard.

Dose Display



## *Radiation Dose Structured Reports*

---

- Radiation Dose Structured Reports (RDSRs) are provided in newer software versions in a defined DICOM format
- They provide the most complete set of information regarding the irradiating events
- The reports are very detailed and require an RDSR viewer for easy visualization of relevant information

Dose Display





## Questions

---

- Please contact the medical physicist providing support for your CT practice, your lead technologist, supervising radiologist or manufacturer's application specialist with questions regarding these important topics and concepts.



## Acknowledgements

---

- **AAPM**
  - Dianna Cody, Dustin Gress, Michael Heard, Jim Kofler, Cynthia McCollough, Mike McNitt-Gray, Bob Pizzutiello, Mark Supanich
- **ACR**
  - Mark Armstrong, Penny Butler, Dina Hernandez
- **ASRT**
  - Virginia Lester
- **DICOM**
  - David Clunie, Kevin O'Donnell
- **FDA**
  - Thalia Mills

A special thank you to Dr. Mark Supanich for his considerable efforts in leading the working group in developing these slides.



## Acknowledgements

---

- GE
  - John Jaeckle
- Hitachi
  - Mark Silverman
- Philips
  - Amar Dhanantwari
- Neusoft
  - Keith Mildemberger
- NeuroLogica
  - Donald Fickett
- Siemens
  - Christianne Liedecker
- Toshiba
  - Kristen Boedecker
- MITA
  - Brian Abraham